

Composition

| Dulbecco's Phosphate Buffered Saline (DPBS) | | CAS number | R1501 Liquid mg/l | R1502 Liquid mg/l | R1503 Liquid mg/l |
|----------------------------------------------------|-----------------------------------------|------------|-------------------|-------------------|-------------------|
| Inorganic Salts | Calcium Chloride Dihydrate | 10035-04-8 | 1330 | - | - |
| | Magnesium Chloride Hexahydrate | 7791-18-6 | 1000 | - | - |
| | Potassium Chloride | 7447-40-7 | 2000 | 2000 | 200 |
| | Potassium Phosphate Monobasic Anhydrous | 7778-77-0 | 2000 | 2000 | 200 |
| | Sodium Chloride | 7647-14-5 | 80000 | 80000 | 8000 |
| | Sodium Phosphate Dibasic Anhydrous | 7558-79-4 | 11500 | 11500 | 1150 |