



Composition

Hank's Balance Salts Solution		CAS number	R1525 Liquid mg/l	R1526 Liquid mg/l	R1527 Liquid mg/l	R1528 Liquid mg/l
Inorganic Salts	Calcium Chloride Dihydrate	10035-04-8	185	185	-	-
	Magnesium Sulfate Anhydrous	7487-88-9	97,67	97,67	-	-
	Potassium Chloride	7447-40-7	400	400	400	400
	Potassium Phosphate Monobasic Anhydrous	7778-77-0	60	60	60	60
	Sodium Bicarbonate	144-55-8	350	350	350	350
	Sodium Chloride	7647-14-5	8000	8000	8000	8000
	Sodium Phosphate Dibasic Anhydrous	7558-79-4	47,88	47,88	47,88	47,88
Other	D-Glucose Anhydrous	50-99-7	1000	1000	1000	1000
	Phenol Red Sodium Salt	34487-61-1	11	-	11	-