

## Composition

<b>Hank's Balance Salts Solution</b>		CAS number	R1530 Liquid mg/l	R1533 Liquid mg/l
Inorganic Salts	Calcium Chloride Dihydrate	10035-04-8	1850	-
	Magnesium Sulfate Anhydrous	7487-88-9	976,7	-
	Potassium Chloride	7447-40-7	4000	4000
	Potassium Phosphate Monobasic Anhydrous	7778-77-0	600	600
	Sodium Bicarbonate	144-55-8	-	-
	Sodium Chloride	7647-14-5	80000	80000
	Sodium Phosphate Dibasic Anhydrous	7558-79-4	478,8	478,8
Other	D-Glucose Anhydrous	50-99-7	10000	10000
	Phenol Red Sodium Salt	34487-61-1	110	-