

Composition

Trypsin		CAS number	R1517 Liquid mg/l	R1518 Liquid mg/l
Inorganic Salts	EDTA Disodium Salt Dihydrate	6381-92-6	200	380
	Potassium Chloride	7447-40-7	400	400
	Potassium Phosphate Monobasic Anhydrous	7778-77-0	60	60
	Sodium Bicarbonate	144-55-8	350	350
	Sodium Chloride	7647-14-5	8000	8000
	Sodium Phosphate Dibasic Anhydrous	7558-79-4	48	48
Other	D-Glucose Anhydrous	50-99-7	1000	1000
	Phenol Red Sodium Salt	34487-61-1	10	10
	Trypsin 1:250	9002-07-7	2500	2500