

Composition

MEM Vitamins	CAS number	M1241 Liquid mg/l
Sodium Chloride	7647-14-5	8500
Choline Chloride	67-48-1	100
D-Ca Pantothenate	137-08-6	100
Folic Acid	59-30-3	100
Myo-Inositol	87-89-8	200
Nicotinamide (Nicotinic acid amide)	98-92-0	100
Pyridoxal Hydrochloride	65-22-5	100
Riboflavin	83-88-5	10
Thiamine Hydrochloride	67-03-8	100