

Composition

Antibiotics		CAS Number	A1001 Liquid mg/l	A1002 Liquid mg/l	A1003 Liquid mg/l	A1004 Liquid mg/l	A1005 Liquid mg/l	A1006 Liquid mg/l	A1007 Liquid mg/l
Antibiotics	Amphotericin B	1397-89-3	250	25	-	-	-	-	-
	Penicillin G Sodium Salt	69-57-8	-	6027	-	-	-	6027	60,27
	Streptomycin Sulfate	3810-74-0	-	10000	-	-	-	10000	100
	Gentamicin Sulfate	1405-41-0	-	-	-	10000	50000	-	-
	G-418	108321-42-2	-	-	50000	-	-	-	-
Amino Acids	L-Glutamine	56-85-9	-	-	-	-	-	29200	-
Inorganic Salts	Sodium Chloride	7647-14-5	-	8500	-	-	-	-	8500
	Sodium Deoxycholate	302-95-4	205	-	-	-	-	-	-