



Fetal Bovine Serum

FBS PRODUCT HIGHLIGHTS

- Triple 0.1µm sterile filtered
- Free storage for up to 1 year
- Free Heat Inactivation
- Gamma Irradiation upon request
- Ships on dry ice, gel packs or refrigerated truck
- Store frozen at -20 to -10°C
- Samples upon request
- Available in 100mL, 500mL and larger volumes

All of our FBS has fully traceable sourcing, collected or imported and treated in agreement with European regulations.

Broad geographic origins: South Africa, Europe, South America, Chile, Uruguay, Central America to best serve our customs needs.

Enhanced quality testing: endotoxin, hemoglobin, mycoplasma, virus and performance assays.

FEATURES

Consistent quality you can rely on: reduce lot-to-lot variability

Always in Stock: no reservations or long waits required

Naturally pure: complete free from synthetic additives

Versatile Performance: ideal for a wide variety of cell lines

APPLICATIONS

Vaccines: low antibody levels, ideal for measles, mumps, rubella production

Cell Culture: rich in growth-promoting factors for robust cell lines

Cloning: high growth factor supports tissue cloning

PRODUCT OPTIONS

Pinnacle Plus FBS – Our Highest Quality FBS

Endotoxin <5 EU/mL

Hemoglobin <25mg/dL

Selected from Pinnacle lots that exhibit superior cell growth achieving >80% cell growth after 6 days on all 4 cell lines (L929, SP2/0-Ag14, HELA and MRC-5).

Chose Pinnacle Plus for the most demanding applications requiring the highest standards with low variations from lot-to-lot. This consistency helps ensure success with your cell lines while eliminates long testing procedures before purchase.

Pinnacle FBS – Premium Quality FBS

Endotoxin <5 EU/mL

Hemoglobin <25mg/dL

Pinnacle FBS has low hemoglobina and endotoxin for reliable and consistent cell growth at a more economical price point. Promotes cell growth on wide range of sensitive cell lines (L929., SP2/0-Ag14, HELA, and MRC-5) with reproducible results.



Phoenix Scientific

For laboratory research and development use only. Not for human use.

PHOENIX SCIENTIFIC

Rua Afonso Almeida nº16, V11

8600-674 Lagos, Portugal

+351 308 800 210

sales-eu@phoenix-scientific.com

www.phoenix-scientific.com